## Suspension freezing as a novel approach for increasing the efficiency of the Laser-Induced Breakdown Spectroscopy method in the study of nano and submicron particles

Shevchenko M.A.<sup>1,@</sup> and Voronova V.<sup>1</sup>

A new method for enhancing Laser-Induced Breakdown Spectroscopy (LIBS) [1] signal intensity in elemental analysis of nano- and sub-micron particles in suspension has been demonstrated. The method is based on the displacement of particles by a solidification front during the process of directed freezing of a liquid [2]. As a result, a densely packed layer of particles is formed on the surface, which leads to an increase in the intensity of plasma emission when interacting with laser radiation. It is shown that the morphology of the resulting layer of particles can be controlled by changing the freezing configuration. The results of increasing the sensitivity of the LIBS method for various suspension freezing configurations, as well as the influence of particle concentration on the efficiency of the process are presented.

<sup>&</sup>lt;sup>1</sup> Lebedev Physical Institute of the Russian Academy of Sciences, Leninsky Avenue 53, Moscow, 119991, Russia

<sup>&</sup>lt;sup>@</sup> mishev87@mail.ru

<sup>[1]</sup> El Haddad J, Canioni L and Bousquet B 2014 Spectrochimica Acta Part B: Atomic Spectroscopy 101 171–182

<sup>[2]</sup> You J, Wang J, Wang L, Wang Z, Wang Z, Li J and Lin X 2017 Colloids and Surfaces A: Physicochemical and Engineering Aspects 531 91–98